

NGFR/TNFRSF16 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51068-1**

Capture Antibody Information

Catalog Number:
60739-1-PBS
Host:
Mouse
Isotype:
IgG2b
Purification Method:
Protein A purification

Clone ID:
1G10A6
Reactivity:
human
Immunogen Catalog Number:
Ag29849

Conjugate:
Unconjugated
Full name:
nerve growth factor receptor (TNFR
superfamily, member 16)
Gene ID:
4804

Detection Antibody Information

Catalog Number:
60739-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2D8A6
Reactivity:
human
GenBank:
NM_002507
Immunogen Catalog Number:
Ag29849

Conjugate:
Unconjugated
Full name:
nerve growth factor receptor (TNFR
superfamily, member 16)
Gene ID:
4804

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP51068-1 targets NGFR/TNFRSF16 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NGFR/TNFRSF16 Monoclonal antibody, PBS Only (Capture) 60739-1-PBS (1G10A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: NGFR/TNFRSF16 Monoclonal antibody, PBS Only (Detector) 60739-2-PBS (2D8A6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

**This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.**

Selected Validation Data



Cytometric bead array standard curve of MP51068-1, p75NTR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60739-1-PBS. Detection antibody: 60739-2-PBS. Standard:Ag29849. Range: 0.781-50 ng/mL