

# INADL Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51174-1**

## Capture Antibody Information

Catalog Number: 67833-2-PBS	Clone ID: 2C9F8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: InaD-like (Drosophila)
Isotype: IgG1	Immunogen Catalog Number: Ag15953	Gene ID: 10207
Purification Method: Protein G Magarose purification		

## Detection Antibody Information

Catalog Number: 67833-3-PBS	Clone ID: 2C7E6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: InaD-like (Drosophila)
Isotype: IgG1	GenBank: BC142661	Gene ID: 10207
Purification Method: Protein G purification	Immunogen Catalog Number: Ag15953	

## Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP51174-1 targets INADL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: INADL Monoclonal antibody, PBS Only (Capture) 67833-2-PBS (2C9F8). 100 µg. Concentration 1 mg/mL.

Detection antibody: INADL Monoclonal antibody, PBS Only (Detector) 67833-3-PBS (2C7E6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51174-1, INADL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67833-2-PBS. Detection antibody: 67833-3-PBS. Standard:Ag15953. Range: 0.195-100 ng/mL