

# MET Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51183-1**

## Capture Antibody Information

**Catalog Number:**  
60813-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2b  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
1B1A10  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag23140

**Conjugate:**  
Unconjugated  
**Full name:**  
met proto-oncogene (hepatocyte growth factor receptor)  
**Gene ID:**  
4233

## Detection Antibody Information

**Catalog Number:**  
60813-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2b  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1A2D9  
**Reactivity:**  
human  
**GenBank:**  
BC130420  
**Immunogen Catalog Number:**  
Ag23140

**Conjugate:**  
Unconjugated  
**Full name:**  
met proto-oncogene (hepatocyte growth factor receptor)  
**Gene ID:**  
4233

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51183-1 targets MET in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MET Monoclonal antibody, PBS Only (Capture) 60813-1-PBS (1B1A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: MET Monoclonal antibody, PBS Only (Detector) 60813-2-PBS (1A2D9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP51183-1, MET Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60813-1-PBS. Detection antibody: 60813-2-PBS. Standard:Ag23140. Range: 0.781-50 ng/mL