

For Research Use Only

# TULP2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51190-2

## Capture Antibody Information

Catalog Number:  
60818-3-PBS  
Host:  
Mouse  
Isotype:  
IgG1  
Purification Method:  
Protein G Magarose purification

Clone ID:  
1B2E7  
Reactivity:  
human  
Immunogen Catalog Number:  
Ag33166

Conjugate:  
Unconjugated  
Full name:  
tubby like protein 2  
Gene ID:  
7288

## Detection Antibody Information

Catalog Number:  
60818-4-PBS  
Host:  
Mouse  
Isotype:  
IgG1  
Purification Method:  
Protein G Magarose purification

Clone ID:  
1F7F1  
Reactivity:  
human  
GenBank:  
BC026070  
Immunogen Catalog Number:  
Ag33166

Conjugate:  
Unconjugated  
Full name:  
tubby like protein 2  
Gene ID:  
7288

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51190-2 targets TULP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TULP2 Monoclonal antibody, PBS Only (Capture) 60818-3-PBS (1B2E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: TULP2 Monoclonal antibody, PBS Only (Detector) 60818-4-PBS (1F7F1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51190-2, TULP2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60818-3-PBS. Detection antibody: 60818-4-PBS. Standard:Ag33166. Range: 0.391-50 ng/mL