

For Research Use Only



IL-10RA Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51203-4

Capture Antibody Information

Catalog Number: 60755-5-PBS	Clone ID: 2F3A8	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 10 receptor, alpha
Isotype: IgG1	Immunogen Catalog Number: Ag4156	Gene ID: 3587
Purification Method: Protein G Magarose purification		

Detection Antibody Information

Catalog Number: 60755-2-PBS	Clone ID: 3B10B4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 10 receptor, alpha
Isotype: IgG1	GenBank: BC028082	Gene ID: 3587
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag4156	

Applications

Tested Applications: Cytometric bead array	Range: 0.098-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51203-4 targets IL-10RA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL-10RA Monoclonal antibody, PBS Only (Capture) 60755-5-PBS (2F3A8). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL-10RA Monoclonal antibody, PBS Only (Detector) 60755-2-PBS (3B10B4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51203-4, IL-10RA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60755-5-PBS. Detection antibody: 60755-2-PBS. Standard:Ag4156. Range: 0.098-100 ng/mL.