

For Research Use Only



RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51247-1

Capture Antibody Information

Catalog Number:
67502-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A5B3

Reactivity:
human

Immunogen Catalog Number:
Ag5782

Conjugate:
Unconjugated

Full name:
Rho-related BTB domain containing 3

Gene ID:
22836

Detection Antibody Information

Catalog Number:
67502-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1G4F9

Reactivity:
human

GenBank:
BC041337

Immunogen Catalog Number:
Ag5782

Conjugate:
Unconjugated

Full name:
Rho-related BTB domain containing 3

Gene ID:
22836

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51247-1 targets RHOBTB3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RHOBTB3 Monoclonal antibody, PBS Only (Capture) 67502-2-PBS (1A5B3). 100 µg. Concentration 1 mg/mL.

Detection antibody: RHOBTB3 Monoclonal antibody, PBS Only (Detector) 67502-3-PBS (1G4F9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51247-1, RHOBTB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67502-2-PBS. Detection antibody: 67502-3-PBS. Standard:Ag5782. Range: 0.391-100 ng/mL