

RECK Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51252-3**

Capture Antibody Information

Catalog Number:
60843-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2F1C5
Reactivity:
human
Immunogen Catalog Number:
Ag26143

Conjugate:
Unconjugated
Full name:
reversion-inducing-cysteine-rich protein with kazal motifs
Gene ID:
8434

Detection Antibody Information

Catalog Number:
60843-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E7G9
Reactivity:
human
GenBank:
BC060806
Immunogen Catalog Number:
Ag26143

Conjugate:
Unconjugated
Full name:
reversion-inducing-cysteine-rich protein with kazal motifs
Gene ID:
8434

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51252-3 targets RECK in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RECK Monoclonal antibody, PBS Only (Capture) 60843-4-PBS (2F1C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: RECK Monoclonal antibody, PBS Only (Detector) 60843-5-PBS (2E7G9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51252-3, RECK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60843-4-PBS. Detection antibody: 60843-5-PBS. Standard:Ag26143. Range: 0.098-12.5 ng/mL