

For Research Use Only



# IL22RA1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51277-3

## Capture Antibody Information

Catalog Number:  
60866-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1F8H3

Reactivity:  
human

Immunogen Catalog Number:  
Ag29983

Conjugate:  
Unconjugated

Full name:  
interleukin 22 receptor, alpha 1

Gene ID:  
58985

## Detection Antibody Information

Catalog Number:  
60866-4-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A Magarose purification

Clone ID:  
2E12B4

Reactivity:  
human

GenBank:  
BC029273

Immunogen Catalog Number:  
Ag29983

Conjugate:  
Unconjugated

Full name:  
interleukin 22 receptor, alpha 1

Gene ID:  
58985

## Applications

Tested Applications:  
Cytometric bead array

Range:  
3.125-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51277-3 targets IL22RA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL22RA1 Monoclonal antibody, PBS Only (Capture/Detector) 60866-2-PBS (1F8H3). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL22RA1 Monoclonal antibody, PBS Only (Detector) 60866-4-PBS (2E12B4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51277-3, IL22RA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60866-2-PBS. Detection antibody: 60866-4-PBS. Standard:Ag29983. Range: 3.125-200 ng/mL.