

GGA1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51292-1**

Capture Antibody Information

Catalog Number:
60864-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C10B4
Reactivity:
human
Immunogen Catalog Number:
Ag22614

Conjugate:
Unconjugated
Full name:
golgi associated, gamma adaptin ear containing, ARF binding protein 1
Gene ID:
26088

Detection Antibody Information

Catalog Number:
60864-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C12C1
Reactivity:
human
GenBank:
BC044629
Immunogen Catalog Number:
Ag22614

Conjugate:
Unconjugated
Full name:
golgi associated, gamma adaptin ear containing, ARF binding protein 1
Gene ID:
26088

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-12.5 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51292-1 targets GGA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GGA1 Monoclonal antibody, PBS Only (Capture) 60864-1-PBS (1C10B4). 100 µg. Concentration 1 mg/mL.

Detection antibody: GGA1 Monoclonal antibody, PBS Only (Detector) 60864-2-PBS (1C12C1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51292-1, GGA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60864-1-PBS. Detection antibody: 60864-2-PBS. Standard:Ag22614. Range: 0.098-12.5 ng/mL