

For Research Use Only



# SPATS2L Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51296-3

## Capture Antibody Information

Catalog Number:  
60881-1-PBS  
Host:  
Mouse  
Isotype:  
IgG1  
Purification Method:  
Protein G purification

Clone ID:  
1C5F4  
Reactivity:  
human  
Immunogen Catalog Number:  
Ag9647

Conjugate:  
Unconjugated  
Full name:  
viral DNA polymerase-transactivated protein 6  
Gene ID:  
26010

## Detection Antibody Information

Catalog Number:  
60881-4-PBS  
Host:  
Mouse  
Isotype:  
IgG2a  
Purification Method:  
Protein A Magarose purification

Clone ID:  
1F6F8  
Reactivity:  
human  
GenBank:  
BC018736  
Immunogen Catalog Number:  
Ag9647

Conjugate:  
Unconjugated  
Full name:  
viral DNA polymerase-transactivated protein 6  
Gene ID:  
26010

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51296-3 targets SPATS2L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPATS2L Monoclonal antibody, PBS Only (Capture) 60881-1-PBS (1C5F4). 100 µg. Concentration 1 mg/mL.

Detection antibody: SPATS2L Monoclonal antibody, PBS Only (Detector) 60881-4-PBS (1F6F8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51296-3, SPATS2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60881-1-PBS. Detection antibody: 60881-4-PBS. Standard:Ag9647. Range: 0.391-100 ng/mL.