

For Research Use Only

BMP-4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51307-1

Capture Antibody Information

Catalog Number:
60897-1-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A purification

Clone ID:
1B8E11

Reactivity:
human

Immunogen Catalog Number:
Ag28840

Conjugate:
Unconjugated

Full name:
bone morphogenetic protein 4

Gene ID:
652

Detection Antibody Information

Catalog Number:
69015-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein A purification

Clone ID:
3B5A7

Reactivity:
human

GenBank:
Immunogen Catalog Number:

HZ-1045

Conjugate:
Unconjugated

Full name:
bone morphogenetic protein 4

Gene ID:
652

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51307-1 targets BMP-4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BMP4 Monoclonal antibody, PBS Only (Capture) 60897-1-PBS (1B8E11). 100 µg. Concentration 1 mg/mL.

Detection antibody: NeutraKine® BMP-4 Monoclonal antibody, PBS Only (Detector) 69015-1-PBS (3B5A7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of cell line supernatant using MP51307-1, BMP-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60897-1-PBS. Detection antibody: 69015-1-PBS. Supernatant (50µL/well) from the following cell lines were assayed: HT-29, TF-1.



Cytometric bead array standard curve of MP51307-1, BMP-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60897-1-PBS. Detection antibody: 69015-1-PBS. Standard:HZ-1045. Range: 6.25-100 ng/mL.