

For Research Use Only



HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51320-1

Capture Antibody Information

Catalog Number:

67445-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

3D2D10

Reactivity:

human

Immunogen Catalog Number:

Ag6699

Conjugate:

Unconjugated

Full name:

heterogeneous nuclear
ribonucleoprotein A2/B1

Gene ID:

3181

Detection Antibody Information

Catalog Number:

67445-3-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1C2B7

Reactivity:

human

GenBank:

BC000506

Immunogen Catalog Number:

Ag6699

Conjugate:

Unconjugated

Full name:

heterogeneous nuclear
ribonucleoprotein A2/B1

Gene ID:

3181

Applications

Tested Applications:

Cytometric bead array

Range:

0.098-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:

It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP51320-1 targets HNRNPA2B1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HNRNPA2B1 Monoclonal antibody, PBS Only (Capture) 67445-2-PBS (3D2D10). 100 µg.
Concentration 1 mg/mL.

Detection antibody: HNRNPA2B1 Monoclonal antibody, PBS Only (Detector) 67445-3-PBS (1C2B7). 100 µg.
Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51320-1, HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67445-2-PBS. Detection antibody: 67445-3-PBS.



Cytometric bead array standard curve of MP51320-1, HNRNPA2B1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67445-2-PBS. Detection antibody: 67445-3-PBS. Standard:Ag6699. Range: 0.098-100 ng/mL.