

# CHAC2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51336-1**

## Capture Antibody Information

**Catalog Number:**  
60914-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
2C11E2  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag9586

**Conjugate:**  
Unconjugated  
**Full name:**  
ChaC, cation transport regulator homolog 2 (E. coli)  
**Gene ID:**  
494143

## Detection Antibody Information

**Catalog Number:**  
60914-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1F4B10  
**Reactivity:**  
human  
**GenBank:**  
BC019239  
**Immunogen Catalog Number:**  
Ag9586

**Conjugate:**  
Unconjugated  
**Full name:**  
ChaC, cation transport regulator homolog 2 (E. coli)  
**Gene ID:**  
494143

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.195-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51336-1 targets CHAC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CHAC2 Monoclonal antibody, PBS Only (Capture) 60914-1-PBS (2C11E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: CHAC2 Monoclonal antibody, PBS Only (Detector) 60914-2-PBS (1F4B10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

Selected Validation Data



Sample test of MP51336-1, CHAC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60914-1-PBS. Detection antibody: 60914-2-PBS.



Cytometric bead array standard curve of MP51336-1, CHAC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60914-1-PBS. Detection antibody: 60914-2-PBS. Standard:Ag9586. Range: 0.195-200 ng/mL