

For Research Use Only

# LCORL Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51368-1

## Capture Antibody Information

**Catalog Number:**  
60954-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
3D5E11  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag11344

**Conjugate:**  
Unconjugated  
**Full name:**  
ligand dependent nuclear receptor corepressor-like  
**Gene ID:**  
254251

## Detection Antibody Information

**Catalog Number:**  
60954-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
3A11F2  
**Reactivity:**  
human  
**GenBank:**  
BC037322  
**Immunogen Catalog Number:**  
Ag11344

**Conjugate:**  
Unconjugated  
**Full name:**  
ligand dependent nuclear receptor corepressor-like  
**Gene ID:**  
254251

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.195-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51368-1 targets LCORL in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LCORL Monoclonal antibody, PBS Only (Capture) 60954-1-PBS (3D5E11). 100 µg. Concentration 1 mg/mL.

Detection antibody: LCORL Monoclonal antibody, PBS Only (Detector) 60954-2-PBS (3A11F2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sample test of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS.

Cytometric bead array standard curve of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS. Standard: Ag11344. Range: 0.195-50 ng/mL.