

For Research Use Only

# TAF11 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51378-1

## Capture Antibody Information

**Catalog Number:**  
60962-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
3B1G7  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag9629

**Conjugate:**  
Unconjugated  
**Full name:**  
TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa  
**Gene ID:**  
6882

## Detection Antibody Information

**Catalog Number:**  
60962-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG2a  
**Purification Method:**  
Protein A Magarose purification

**Clone ID:**  
3B10G9  
**Reactivity:**  
human  
**GenBank:**  
BC021972  
**Immunogen Catalog Number:**  
Ag9629

**Conjugate:**  
Unconjugated  
**Full name:**  
TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa  
**Gene ID:**  
6882

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
12.5-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51378-1 targets TAF11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAF11 Monoclonal antibody, PBS Only (Capture) 60962-1-PBS (3B1G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: TAF11 Monoclonal antibody, PBS Only (Detector) 60962-2-PBS (3B10G9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.



## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sample test of MP51378-1, TAF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60962-1-PBS. Detection antibody: 60962-2-PBS.

Cytometric bead array standard curve of MP51378-1, TAF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60962-1-PBS. Detection antibody: 60962-2-PBS. Standard: Ag9629. Range: 12.5-200 ng/mL.