

DDHD2 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51382-2**

Capture Antibody Information

Catalog Number:
60965-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1D9A1
Reactivity:
human
Immunogen Catalog Number:
Ag18391

Conjugate:
Unconjugated
Full name:
DDHD domain containing 2
Gene ID:
23259

Detection Antibody Information

Catalog Number:
60965-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3E1C9
Reactivity:
human
GenBank:
BC010504
Immunogen Catalog Number:
Ag18391

Conjugate:
Unconjugated
Full name:
DDHD domain containing 2
Gene ID:
23259

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51382-2 targets DDHD2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDHD2 Monoclonal antibody, PBS Only (Capture) 60965-3-PBS (1D9A1). 100 µg. Concentration 1 mg/mL.

Detection antibody: DDHD2 Monoclonal antibody, PBS Only (Detector) 60965-2-PBS (3E1C9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51382-2, DDHD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60965-3-PBS. Detection antibody: 60965-2-PBS. Standard:Ag18391. Range: 0.391-100 ng/mL.