

For Research Use Only

RALY Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51385-1

Capture Antibody Information

Catalog Number:

68011-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

3A10D7

Reactivity:

human

Immunogen Catalog Number:

Ag12250

Conjugate:

Unconjugated

Full name:

RNA binding protein, autoantigenic (hnRNP-associated with lethal yellow homolog (mouse))

Gene ID:

22913

Detection Antibody Information

Catalog Number:

68011-3-PBS

Host:

Mouse

Isotype:

IgG3

Purification Method:

Protein A Magarose purification

Clone ID:

2H1H10

Reactivity:

human

GenBank:

BC103753

Immunogen Catalog Number:

Ag12250

Conjugate:

Unconjugated

Full name:

RNA binding protein, autoantigenic (hnRNP-associated with lethal yellow homolog (mouse))

Gene ID:

22913

Applications

Tested Applications:

Cytometric bead array

Range:

0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51385-1 targets RALY in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RALY Monoclonal antibody, PBS Only (Capture) 68011-2-PBS (3A10D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: RALY Monoclonal antibody, PBS Only (Detector) 68011-3-PBS (2H1H10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

