

p65; RELA Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51387-1**

Capture Antibody Information

Catalog Number:
66535-2-PBS
Host:
Mouse
Isotype:
IgM
Purification Method:
Euglobulin precipitation

Clone ID:
2B6G11
Reactivity:
human

Conjugate:
Unconjugated
Full name:
v-rel reticuloendotheliosis viral oncogene homolog A (avian)
Gene ID:
5970

Detection Antibody Information

Catalog Number:
66535-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1D8H3
Reactivity:
human
GenBank:
BC011603
Immunogen Catalog Number:
Ag1199

Conjugate:
Unconjugated
Full name:
v-rel reticuloendotheliosis viral oncogene homolog A (avian)
Gene ID:
5970

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51387-1 targets p65; RELA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: p65; RELA Monoclonal antibody, PBS Only (Capture) 66535-2-PBS (2B6G11). 100 µg. Concentration 1 mg/mL.

Detection antibody: p65; RELA Monoclonal antibody, PBS Only (Detector) 66535-3-PBS (1D8H3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Sample test of MP51387-1, p65; RELA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66535-2-PBS. Detection antibody: 66535-3-PBS.



Cytometric bead array standard curve of MP51387-1, p65; RELA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66535-2-PBS. Detection antibody: 66535-3-PBS. Standard: Ag1199. Range: 1.563-200 ng/mL.