

For Research Use Only

TAF13 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51418-1

Capture Antibody Information

Catalog Number:
60984-1-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
2G4C12
Reactivity:
human
Immunogen Catalog Number:
Ag18502

Conjugate:
Unconjugated
Full name:
TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa
Gene ID:
6884

Detection Antibody Information

Catalog Number:
60984-2-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
2A1D8
Reactivity:
human
GenBank:
BC121180
Immunogen Catalog Number:
Ag18502

Conjugate:
Unconjugated
Full name:
TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa
Gene ID:
6884

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51418-1 targets TAF13 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: TAF13 Monoclonal antibody, PBS Only (Capture) 60984-1-PBS (2G4C12). 100 µg. Concentration 1 mg/mL.
Detection antibody: TAF13 Monoclonal antibody, PBS Only (Detector) 60984-2-PBS (2A1D8). 100 µg. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.



Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51418-1, TAF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60984-1-PBS. Detection antibody: 60984-2-PBS. Standard: Ag18502. Range: 0.098-100 ng/mL.

Cytometric bead array sample test of MP51418-1, TAF13 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60984-1-PBS. Detection antibody: 60984-2-PBS.