

For Research Use Only

RNF19A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51425-1

Capture Antibody Information

Catalog Number:
60990-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B8B5

Reactivity:
human

Immunogen Catalog Number:
Ag33826

Conjugate:
Unconjugated

Full name:
ring finger protein 19A

Gene ID:
25897

Detection Antibody Information

Catalog Number:
60990-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2H7F1

Reactivity:
human

GenBank:
BC093938

Immunogen Catalog Number:
Ag33826

Conjugate:
Unconjugated

Full name:
ring finger protein 19A

Gene ID:
25897

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51425-1 targets RNF19A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RNF19A Monoclonal antibody, PBS Only (Capture) 60990-1-PBS (1B8B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: RNF19A Monoclonal antibody, PBS Only (Detector) 60990-2-PBS (2H7F1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51425-1, RNF19A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60990-1-PBS. Detection antibody: 60990-2-PBS. Standard: Ag33826. Range: 0.391-400 ng/mL.



Cytometric bead array sample test of MP51425-1, RNF19A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60990-1-PBS. Detection antibody: 60990-2-PBS.