

For Research Use Only

# NMNAT1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51431-1

## Capture Antibody Information

**Catalog Number:**  
60994-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1B4G9  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag29593

**Conjugate:**  
Unconjugated  
**Full name:**  
nicotinamide nucleotide  
adenyltransferase 1  
**Gene ID:**  
64802

## Detection Antibody Information

**Catalog Number:**  
60994-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG3  
**Purification Method:**  
Protein A purification

**Clone ID:**  
1G7B7  
**Reactivity:**  
human  
**GenBank:**  
BC014943  
**Immunogen Catalog Number:**  
Ag29593

**Conjugate:**  
Unconjugated  
**Full name:**  
nicotinamide nucleotide  
adenyltransferase 1  
**Gene ID:**  
64802

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.391-100 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51431-1 targets NMNAT1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: NMNAT1 Monoclonal antibody, PBS Only (Capture) 60994-1-PBS (1B4G9). 100 µg. Concentration 1 mg/mL.  
Detection antibody: NMNAT1 Monoclonal antibody, PBS Only (Detector) 60994-2-PBS (1G7B7). 100 µg. Concentration 1 mg/mL.  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.



## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51431-1, NMNAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60994-1-PBS. Detection antibody: 60994-2-PBS. Standard: Ag29593. Range: 0.391-100 ng/mL.

Cytometric bead array sample test of MP51431-1, NMNAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60994-1-PBS. Detection antibody: 60994-2-PBS.