

For Research Use Only

IRF5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51481-1

Capture Antibody Information

Catalog Number:
66835-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2E1B1
Reactivity:
human
Immunogen Catalog Number:
Ag5442

Conjugate:
Unconjugated
Full name:
interferon regulatory factor 5
Gene ID:
3663

Detection Antibody Information

Catalog Number:
66835-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3D1E11
Reactivity:
human
GenBank:
BC004201
Immunogen Catalog Number:
Ag5442

Conjugate:
Unconjugated
Full name:
interferon regulatory factor 5
Gene ID:
3663

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51481-1 targets IRF5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IRF5 Monoclonal antibody, PBS Only (Capture) 66835-2-PBS (2E1B1). 100 µg. Concentration 1 mg/mL.

Detection antibody: IRF5 Monoclonal antibody, PBS Only (Detector) 66835-3-PBS (3D1E11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

