

For Research Use Only

# RBL1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51495-1

## Capture Antibody Information

Catalog Number:  
61017-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1C7D11

Reactivity:  
human

Immunogen Catalog Number:  
Ag23612

Conjugate:  
Unconjugated

Full name:  
retinoblastoma-like 1 (p107)

Gene ID:  
5933

## Detection Antibody Information

Catalog Number:  
61017-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2D8D5

Reactivity:  
human

GenBank:  
BC032247

Immunogen Catalog Number:  
Ag23612

Conjugate:  
Unconjugated

Full name:  
retinoblastoma-like 1 (p107)

Gene ID:  
5933

## Applications

Tested Applications:  
Cytometric bead array

Range:  
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51495-1 targets RBL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBL1 Monoclonal antibody, PBS Only (Capture) 61017-1-PBS (1C7D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBL1 Monoclonal antibody, PBS Only (Detector) 61017-2-PBS (2D8D5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51495-1, RBL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61017-1-PBS. Detection antibody: 61017-2-PBS.



Cytometric bead array standard curve of MP51495-1, RBL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61017-1-PBS. Detection antibody: 61017-2-PBS. Standard: Ag23612. Range: 3.125-400 ng/mL.