

NR3C1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51534-1**

Capture Antibody Information

Catalog Number:
66904-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2E8B5
Reactivity:
human
Immunogen Catalog Number:
Ag28431

Conjugate:
Unconjugated
Full name:
nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)
Gene ID:
2908

Detection Antibody Information

Catalog Number:
66904-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2B10B9
Reactivity:
human
GenBank:
BC015610
Immunogen Catalog Number:
Ag28431

Conjugate:
Unconjugated
Full name:
nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)
Gene ID:
2908

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51534-1 targets NR3C1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Glucocorticoid receptor Monoclonal antibody, PBS Only (Capture) 66904-2-PBS (2E8B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: Glucocorticoid receptor Monoclonal antibody, PBS Only (Detector) 66904-3-PBS (2B10B9). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS.



Cytometric bead array standard curve of MP51534-1, NR3C1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66904-2-PBS. Detection antibody: 66904-3-PBS. Standard:Ag28431. Range: 0.781-100 ng/mL.