

For Research Use Only

FBRS Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51558-1

Capture Antibody Information

Catalog Number:

61040-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

5H5D2

Reactivity:

human

Immunogen Catalog Number:

Ag1870

Conjugate:

Unconjugated

Full name:

fibrosin

Gene ID:

64319

Detection Antibody Information

Catalog Number:

61041-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1C1E1

Reactivity:

human

GenBank:

BC014212

Immunogen Catalog Number:

Ag30051

Conjugate:

Unconjugated

Full name:

fibrosin

Gene ID:

64319

Applications

Tested Applications:

Cytometric bead array

Range:

0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51558-1 targets FBRS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FBRS Monoclonal antibody, PBS Only (Capture) 61040-1-PBS (5H5D2). 100 µg. Concentration 1 mg/mL.

Detection antibody: FBRS Monoclonal antibody, PBS Only (Detector) 61041-1-PBS (1C1E1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51558-1, FBR5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61040-1-PBS. Detection antibody: 61041-1-PBS.



Cytometric bead array standard curve of MP51558-1, FBR5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61040-1-PBS. Detection antibody: 61041-1-PBS. Standard: Ag1870. Range: 0.781-100 ng/mL