

For Research Use Only

# IKZF1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51595-1

## Capture Antibody Information

Catalog Number:  
66966-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1C1G12

Reactivity:  
human

Immunogen Catalog Number:  
Ag28637

Conjugate:  
Unconjugated

Full name:  
IKAROS family zinc finger 1 (Ikaros)

Gene ID:  
10320

## Detection Antibody Information

Catalog Number:  
66966-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
2C9G5

Reactivity:  
human

GenBank:  
BC018349

Immunogen Catalog Number:  
Ag28637

Conjugate:  
Unconjugated

Full name:  
IKAROS family zinc finger 1 (Ikaros)

Gene ID:  
10320

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51595-1 targets IKZF1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IKZF1 Monoclonal antibody, PBS Only (Capture) 66966-2-PBS (1C1G12). 100 µg. Concentration 1 mg/mL.

Detection antibody: IKZF1 Monoclonal antibody, PBS Only (Detector) 66966-3-PBS (2C9G5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS.



Cytometric bead array standard curve of MP51595-1, IKZF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66966-2-PBS. Detection antibody: 66966-3-PBS. Standard: Ag28637. Range: 1.563-200 ng/mL.