

For Research Use Only

SMURF2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51609-1

Capture Antibody Information

Catalog Number:
61048-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1B12B2

Reactivity:
human

Immunogen Catalog Number:
Ag25734

Conjugate:
Unconjugated

Full name:
SMAD specific E3 ubiquitin protein
ligase 2

Gene ID:
64750

Detection Antibody Information

Catalog Number:
61048-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1F6D1

Reactivity:
human

GenBank:
BC093876

Immunogen Catalog Number:
Ag25734

Conjugate:
Unconjugated

Full name:
SMAD specific E3 ubiquitin protein
ligase 2

Gene ID:
64750

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP51609-1 targets SMURF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMURF2 Monoclonal antibody, PBS Only (Capture) 61048-1-PBS (1B12B2). 100 µg. Concentration 1 mg/mL.

Detection antibody: SMURF2 Monoclonal antibody, PBS Only (Detector) 61048-2-PBS (1F6D1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51609-1, SMURF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61048-1-PBS. Detection antibody: 61048-2-PBS.



Cytometric bead array standard curve of MP51609-1, SMURF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61048-1-PBS. Detection antibody: 61048-2-PBS. Standard: Ag25734. Range: 0.781-100 ng/mL