

For Research Use Only

TMSB4X Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51621-1

Capture Antibody Information

Catalog Number:

61052-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1F1C4

Reactivity:

human

Immunogen Catalog Number:

Ag13914

Conjugate:

Unconjugated

Full name:

thymosin beta 4, X-linked

Gene ID:

7114

Detection Antibody Information

Catalog Number:

61053-1-PBS

Host:

Mouse

Isotype:

IgG3

Purification Method:

Protein A Magarose purification

Clone ID:

3D3G1

Reactivity:

human

GenBank:

BC001631

Immunogen Catalog Number:

Ag14097

Conjugate:

Unconjugated

Full name:

thymosin beta 4, X-linked

Gene ID:

7114

Applications

Tested Applications:

Cytometric bead array

Range:

0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51621-1 targets TMSB4X in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMSB4X Monoclonal antibody, PBS Only (Capture) 61052-1-PBS (1F1C4). 100 µg. Concentration 1 mg/mL.

Detection antibody: TMSB4X Monoclonal antibody, PBS Only (Detector) 61053-1-PBS (3D3G1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

