

For Research Use Only

# PFKFB1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51639-1

## Capture Antibody Information

Catalog Number:

61064-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

1A5E2

Reactivity:

human

Immunogen Catalog Number:

Ag16752

Conjugate:

Unconjugated

Full name:

6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 1

Gene ID:

5207

## Detection Antibody Information

Catalog Number:

61064-2-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1H7E3

Reactivity:

human

GenBank:

BC096079

Immunogen Catalog Number:

Ag16752

Conjugate:

Unconjugated

Full name:

6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 1

Gene ID:

5207

## Applications

Tested Applications:

Cytometric bead array

Range:

3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51639-1 targets PFKFB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PFKFB1 Monoclonal antibody, PBS Only (Capture) 61064-1-PBS (1A5E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: PFKFB1 Monoclonal antibody, PBS Only (Detector) 61064-2-PBS (1H7E3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)


W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51639-1, PFKFB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61064-1-PBS. Detection antibody: 61064-2-PBS.



Cytometric bead array standard curve of MP51639-1, PFKFB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61064-1-PBS. Detection antibody: 61064-2-PBS. Standard: Ag16752. Range: 3.125-400 ng/mL