

For Research Use Only

TPM4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51644-1

Capture Antibody Information

Catalog Number:
67244-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1E11G8
Reactivity:
human
Immunogen Catalog Number:
Ag6947

Conjugate:
Unconjugated
Full name:
tropomyosin 4
Gene ID:
7171

Detection Antibody Information

Catalog Number:
67244-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2E12D12
Reactivity:
human, mouse, rat, pig, rabbit
GenBank:
BC037576
Immunogen Catalog Number:
Ag6947

Conjugate:
Unconjugated
Full name:
tropomyosin 4
Gene ID:
7171

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51644-1 targets TPM4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: TPM4 Monoclonal antibody, PBS Only (Capture) 67244-2-PBS (1E11G8). 100 µg. Concentration 1 mg/mL.
Detection antibody: TPM4 Monoclonal antibody, PBS Only (Detector) 67244-1-PBS (2E12D12). 100 µg. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.



Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51644-1, TPM4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67244-2-PBS. Detection antibody: 67244-1-PBS.

Cytometric bead array standard curve of MP51644-1, TPM4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67244-2-PBS. Detection antibody: 67244-1-PBS. Standard: Ag6947. Range: 1.563-200 ng/mL.