

For Research Use Only

BSEP Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51654-1

Capture Antibody Information

Catalog Number:
67512-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2C7C1
Reactivity:
human
Immunogen Catalog Number:
Ag29135

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family B (MDR/TAP), member 11
Gene ID:
8647

Detection Antibody Information

Catalog Number:
67512-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3B3E7
Reactivity:
human
GenBank:
NM_003742
Immunogen Catalog Number:
Ag29135

Conjugate:
Unconjugated
Full name:
ATP-binding cassette, sub-family B (MDR/TAP), member 11
Gene ID:
8647

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51654-1 targets BSEP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BSEP Monoclonal antibody, PBS Only (Capture) 67512-2-PBS (2C7C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: BSEP Monoclonal antibody, PBS Only (Detector) 67512-3-PBS (3B3E7). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

