

For Research Use Only

TAF6L Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51658-1

Capture Antibody Information

Catalog Number:
67569-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2F3H9
Reactivity:
human
Immunogen Catalog Number:
Ag8202

Conjugate:
Unconjugated
Full name:
TAF6-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa
Gene ID:
10629

Detection Antibody Information

Catalog Number:
67569-3-PBS
Host:
Mouse
Isotype:
IgG2a
Purification Method:
Protein A Magarose purification

Clone ID:
1F10E9
Reactivity:
human
GenBank:
BC008785
Immunogen Catalog Number:
Ag8202

Conjugate:
Unconjugated
Full name:
TAF6-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa
Gene ID:
10629

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51658-1 targets TAF6L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: TAF6L Monoclonal antibody, PBS Only (Capture) 67569-2-PBS (2F3H9). 100 µg. Concentration 1 mg/mL.
Detection antibody: TAF6L Monoclonal antibody, PBS Only (Detector) 67569-3-PBS (1F10E9). 100 µg. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51658-1, TAF6L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67569-2-PBS. Detection antibody: 67569-3-PBS.



Cytometric bead array standard curve of MP51658-1, TAF6L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67569-2-PBS. Detection antibody: 67569-3-PBS. Standard: Ag8180. Range: 1.563-200 ng/mL.