

For Research Use Only

EEF1D Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51681-1

Capture Antibody Information

Catalog Number:
60085-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
3E4C2
Reactivity:
human
Immunogen Catalog Number:
Ag0983

Conjugate:
Unconjugated
Full name:
eukaryotic translation elongation factor 1 delta (guanine nucleotide exchange protein)
Gene ID:
1936

Detection Antibody Information

Catalog Number:
60085-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3B1B11
Reactivity:
human, mouse, rat
GenBank:
BC007847
Immunogen Catalog Number:
Ag0983

Conjugate:
Unconjugated
Full name:
eukaryotic translation elongation factor 1 delta (guanine nucleotide exchange protein)
Gene ID:
1936

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51681-1 targets EEF1D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: EEF1D Monoclonal antibody, PBS Only (Capture) 60085-3-PBS (3E4C2). 100 µg. Concentration 1 mg/mL.
Detection antibody: EEF1D Monoclonal antibody, PBS Only (Detector) 60085-1-PBS (3B1B11). 100 µg. Concentration 1 mg/mL.
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

