

FBXW7 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51687-1**

Capture Antibody Information

Catalog Number:
67234-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1F10C6
Reactivity:
human
Immunogen Catalog Number:
Ag28471

Conjugate:
Unconjugated
Full name:
F-box and WD repeat domain containing 7
Gene ID:
55294

Detection Antibody Information

Catalog Number:
67234-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein A purification

Clone ID:
2A10E5
Reactivity:
human
GenBank:
NM_001349798
Immunogen Catalog Number:
Ag28471

Conjugate:
Unconjugated
Full name:
F-box and WD repeat domain containing 7
Gene ID:
55294

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51687-1 targets FBXW7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FBXW7 Monoclonal antibody, PBS Only (Capture) 67234-2-PBS (1F10C6). 100 µg. Concentration 1 mg/mL.

Detection antibody: FBXW7 Monoclonal antibody, PBS Only (Detector) 67234-1-PBS (2A10E5). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51687-1, FBXW7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67234-2-PBS. Detection antibody: 67234-1-PBS.

