

For Research Use Only

# SLC36A3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51713-1

## Capture Antibody Information

**Catalog Number:**  
67929-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
2C1D1  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag16556

**Conjugate:**  
Unconjugated  
**Full name:**  
solute carrier family 36  
(proton/amino acid symporter),  
member 3  
**Gene ID:**  
285641

## Detection Antibody Information

**Catalog Number:**  
67929-3-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G purification

**Clone ID:**  
1H9D7  
**Reactivity:**  
human  
**GenBank:**  
BC101092  
**Immunogen Catalog Number:**  
Ag16556

**Conjugate:**  
Unconjugated  
**Full name:**  
solute carrier family 36  
(proton/amino acid symporter),  
member 3  
**Gene ID:**  
285641

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.563-200 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51713-1 targets SLC36A3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: SLC36A3 Monoclonal antibody, PBS Only (Capture) 67929-2-PBS (2C1D1). 100 µg. Concentration 1 mg/mL.  
Detection antibody: SLC36A3 Monoclonal antibody, PBS Only (Detector) 67929-3-PBS (1H9D7). 100 µg. Concentration 1 mg/mL.  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

