

For Research Use Only

# CTAK1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51716-1

## Capture Antibody Information

**Catalog Number:**  
61082-1-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1D9H11  
**Reactivity:**  
human  
**Immunogen Catalog Number:**  
Ag8303

**Conjugate:**  
Unconjugated  
**Full name:**  
MAP/microtubule affinity-regulating kinase 3  
**Gene ID:**  
4140

## Detection Antibody Information

**Catalog Number:**  
61082-2-PBS  
**Host:**  
Mouse  
**Isotype:**  
IgG1  
**Purification Method:**  
Protein G Magarose purification

**Clone ID:**  
1E3D3  
**Reactivity:**  
human  
**GenBank:**  
BC024773  
**Immunogen Catalog Number:**  
Ag8303

**Conjugate:**  
Unconjugated  
**Full name:**  
MAP/microtubule affinity-regulating kinase 3  
**Gene ID:**  
4140

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
3.125-400 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51716-1 targets CTAK1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: CTAK1 Monoclonal antibody, PBS Only (Capture) 61082-1-PBS (1D9H11). 100 µg. Concentration 1 mg/mL.  
Detection antibody: CTAK1 Monoclonal antibody, PBS Only (Detector) 61082-2-PBS (1E3D3). 100 µg. Concentration 1 mg/mL.  
Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.  
Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.  
Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51716-1, CTAK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61082-1-PBS. Detection antibody: 61082-2-PBS.



Cytometric bead array standard curve of MP51716-1, CTAK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61082-1-PBS. Detection antibody: 61082-2-PBS. Standard: Ag8303. Range: 3.125-400 ng/mL