

For Research Use Only

# C20orf11 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51718-1

## Capture Antibody Information

Catalog Number:  
61088-1-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
2G12E7

Reactivity:  
human

Immunogen Catalog Number:  
Ag19740

Conjugate:  
Unconjugated

Full name:  
chromosome 20 open reading frame  
11

Gene ID:  
54994

## Detection Antibody Information

Catalog Number:  
61088-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1C6C8

Reactivity:  
human

GenBank:  
BC032120

Immunogen Catalog Number:  
Ag19740

Conjugate:  
Unconjugated

Full name:  
chromosome 20 open reading frame  
11

Gene ID:  
54994

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.156-20 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP51718-1 targets C20orf11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C20orf11 Monoclonal antibody, PBS Only (Capture) 61088-1-PBS (2G12E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: C20orf11 Monoclonal antibody, PBS Only (Detector) 61088-2-PBS (1C6C8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51718-1, C20orf11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61088-1-PBS. Detection antibody: 61088-2-PBS.



Cytometric bead array standard curve of MP51718-1, C20orf11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61088-1-PBS. Detection antibody: 61088-2-PBS. Standard: Ag19740. Range: 0.156-20 ng/mL.