

PAIP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51721-2**

Capture Antibody Information

Catalog Number:
68042-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G4B4
Reactivity:
human
Immunogen Catalog Number:
Ag16773

Conjugate:
Unconjugated
Full name:
poly(A) binding protein interacting protein 1
Gene ID:
10605

Detection Antibody Information

Catalog Number:
68042-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1C2F12
Reactivity:
human
GenBank:
BC005295
Immunogen Catalog Number:
Ag16773

Conjugate:
Unconjugated
Full name:
poly(A) binding protein interacting protein 1
Gene ID:
10605

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51721-2 targets PAIP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PAIP1 Monoclonal antibody, PBS Only (Capture) 68042-4-PBS (1G4B4). 100 µg. Concentration 1 mg/mL.

Detection antibody: PAIP1 Monoclonal antibody, PBS Only (Detector) 68042-5-PBS (1C2F12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51721-2, PAIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68042-4-PBS. Detection antibody: 68042-5-PBS.



Cytometric bead array standard curve of MP51721-2, PAIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68042-4-PBS. Detection antibody: 68042-5-PBS. Standard: Ag16773. Range: 0.781-50 ng/mL.