

For Research Use Only

PXR Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51736-1

Capture Antibody Information

Catalog Number:
67912-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein A purification

Clone ID:
1G3A7

Reactivity:
human, mouse, rat

Immunogen Catalog Number:
Ag7991

Conjugate:
Unconjugated

Full name:
nuclear receptor subfamily 1, group 1, member 2

Gene ID:
8856

Detection Antibody Information

Catalog Number:
67912-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
1G7D8

Reactivity:
human

GenBank:
BC017304

Immunogen Catalog Number:
Ag7991

Conjugate:
Unconjugated

Full name:
nuclear receptor subfamily 1, group 1, member 2

Gene ID:
8856

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-400 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51736-1 targets PXR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PXR Monoclonal antibody, PBS Only (Capture) 67912-1-PBS (1G3A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: PXR Monoclonal antibody, PBS Only (Detector) 67912-2-PBS (1G7D8). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

