

For Research Use Only



CAMK2D Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51748-1

Capture Antibody Information

Catalog Number:
61090-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G1D4
Reactivity:
human
Immunogen Catalog Number:
Ag7323

Conjugate:
Unconjugated
Full name:
calcium/calmodulin-dependent
protein kinase II delta
Gene ID:
817

Detection Antibody Information

Catalog Number:
61090-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2H5H1
Reactivity:
human
GenBank:
BC032784
Immunogen Catalog Number:
Ag7323

Conjugate:
Unconjugated
Full name:
calcium/calmodulin-dependent
protein kinase II delta
Gene ID:
817

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51748-1 targets CAMK2D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CAMK2D Monoclonal antibody, PBS Only (Capture) 61090-1-PBS (1G1D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: CAMK2D Monoclonal antibody, PBS Only (Detector) 61090-2-PBS (2H5H1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51748-1, CAMK2D Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61090-1-PBS. Detection antibody: 61090-2-PBS.



Cytometric bead array standard curve of MP51748-1, CAMK2D Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61090-1-PBS. Detection antibody: 61090-2-PBS. Standard: Ag7323. Range: 0.781-100 ng/mL.