

C6orf108 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51773-1**

Capture Antibody Information

Catalog Number:
68310-3-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2D5G2
Reactivity:
human
Immunogen Catalog Number:
Ag22382

Conjugate:
Unconjugated
Full name:
chromosome 6 open reading frame 108
Gene ID:
10591

Detection Antibody Information

Catalog Number:
68310-4-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1D5C3
Reactivity:
human
GenBank:
BC011683
Immunogen Catalog Number:
Ag22382

Conjugate:
Unconjugated
Full name:
chromosome 6 open reading frame 108
Gene ID:
10591

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51773-1 targets C6orf108 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C6orf108 Monoclonal antibody, PBS Only (Capture) 68310-3-PBS (2D5G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: C6orf108 Monoclonal antibody, PBS Only (Detector) 68310-4-PBS (1D5C3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51773-1, C6orf108 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68310-3-PBS. Detection antibody: 68310-4-PBS.



Cytometric bead array standard curve of MP51773-1, C6orf108 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68310-3-PBS. Detection antibody: 68310-4-PBS. Standard:Ag22382. Range: 0.781-100 ng/mL