

EIF3H Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51804-1**

Capture Antibody Information

Catalog Number:
61110-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
1B9G3
Reactivity:
human
Immunogen Catalog Number:
Ag30831

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 3, subunit H
Gene ID:
8667

Detection Antibody Information

Catalog Number:
61110-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1G9D1
Reactivity:
human
GenBank:
BC000386
Immunogen Catalog Number:
Ag30831

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 3, subunit H
Gene ID:
8667

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51804-1 targets EIF3H in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EIF3H Monoclonal antibody, PBS Only (Capture) 61110-1-PBS (1B9G3). 100 µg. Concentration 1 mg/mL.

Detection antibody: EIF3H Monoclonal antibody, PBS Only (Detector) 61110-2-PBS (1G9D1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51804-1, EIF3H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61110-1-PBS. Detection antibody: 61110-2-PBS.



Cytometric bead array standard curve of MP51804-1, EIF3H Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61110-1-PBS. Detection antibody: 61110-2-PBS. Standard:Ag30831. Range: 0.781-100 ng/mL