

SNCB Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51827-2**

Capture Antibody Information

Catalog Number:
68046-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2B10E10
Reactivity:
human, mouse, rat, pig, rabbit, chicken
Immunogen Catalog Number:
Ag0772

Conjugate:
Unconjugated
Full name:
synuclein, beta
Gene ID:
6620

Detection Antibody Information

Catalog Number:
68046-3-PBS
Host:
Mouse
Isotype:
IgG3
Purification Method:
Protein A purification

Clone ID:
2A9E3
Reactivity:
human
GenBank:
BC002902
Immunogen Catalog Number:
Ag0772

Conjugate:
Unconjugated
Full name:
synuclein, beta
Gene ID:
6620

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51827-2 targets SNCB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SNCB Monoclonal antibody, PBS Only (Capture) 68046-1-PBS (2B10E10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SNCB Monoclonal antibody, PBS Only (Detector) 68046-3-PBS (2A9E3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51827-2, SNCB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68046-1-PBS. Detection antibody: 68046-3-PBS.



Cytometric bead array standard curve of MP51827-2, SNCB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68046-1-PBS. Detection antibody: 68046-3-PBS. Standard: Ag0772. Range: 0.781-100 ng/mL.