

SPDEF Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51889-1**

Capture Antibody Information

Catalog Number:
68716-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2B4E9
Reactivity:
human
Immunogen Catalog Number:
Ag33739

Conjugate:
Unconjugated
Full name:
SAM pointed domain containing ets transcription factor
Gene ID:
25803

Detection Antibody Information

Catalog Number:
68716-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2A1C1
Reactivity:
human
GenBank:
BC021299
Immunogen Catalog Number:
Ag33739

Conjugate:
Unconjugated
Full name:
SAM pointed domain containing ets transcription factor
Gene ID:
25803

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51889-1 targets SPDEF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SPDEF Monoclonal antibody, PBS Only (Capture) 68716-2-PBS (2B4E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: SPDEF Monoclonal antibody, PBS Only (Detector) 68716-1-PBS (2A1C1). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array sample test of MP51889-1, SPDEF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68716-2-PBS. Detection antibody: 68716-1-PBS.



Cytometric bead array standard curve of MP51889-1, SPDEF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68716-2-PBS. Detection antibody: 68716-1-PBS. Standard:Ag33739. Range: 0.781-100 ng/mL.