

Speedy™ Mouse/Rat HGF One-Step ELISA Kit Datasheet

Please read it entirely before use

Catalogue Number: SE60014

Size: 96T

Sensitivity: 7.8 pg/mL

Range: 62.5-4000 pg/mL

Usage: For the quantitative detection of mouse/rat HGF concentrations in mouse serum, mouse plasma, rat serum, rat plasma, cell culture supernatant and mouse tissue homogenate.


This product is for research use only and not for use in human or animal therapeutic or diagnostic.

Table of content	page
1. Background	3
2. Principle	3
3. Required Materials	3
4. Kit Components and Storage	4
5. Safety Notes	4
6. Sample Collection and Storage	4
7. Regent Preparation	5
8. Assay Procedure Summary	6
9. Validation Data	7
9.1 Standard curve	7
9.2 Precision	8
9.3 Recovery	8
9.4 Sample values	9
9.5 Sensitivity	9
9.6 Linearity	10
9.7 Specificity	10
10. References	10

1. Background

HGF (hepatocyte growth factor) is a potent mitogen for hepatocytes, also known as hepatopoietin A (HBGA). HGF and its receptor MET are involved in cell proliferation, movement, differentiation, invasion, angiogenesis, and apoptosis by activating multiple downstream signaling pathways. HGF/cMet also can serve as biomarkers for gastrointestinal cancers. HGF regulates inflammation including collagen-induced arthritis, autoimmune neuroinflammation, inflammatory bowel disease and airway inflammation. Elevated levels of HGF of tumors and plasma were observed such as in invasive breast carcinoma, glioma, multiple myeloma and sarcomas.

2. Principle



An anti-tag antibody is pre-coated onto the bottom of wells. After adding antigen or samples, Tagged-Capture antibody and HRP-conjugated detection antibody, a sandwich complex is formed in the solution. TMB acts as a HRP substrate, and the solution color will change from colorless to blue. A stop solution containing sulfuric acid turns the solution yellow. The color intensity is proportional to the quantity of bound protein, which is measurable at 450 nm with the correction wavelength set at 630 nm.

3. Required Materials

- 3.1 A microplate reader capable of measuring absorbance at 450 nm with the correction wavelength set at 630 nm.
- 3.2 Calibrated, adjustable precision pipettes and disposable plastic tips. A manifold multi-channel pipette is recommended for large assays.
- 3.3 Plate washer: automated or manual.
- 3.4 Absorbent paper towels.
- 3.5 Glass or plastic tubes to prepare standard and sample dilutions.
- 3.6 Beakers and graduated cylinders.
- 3.7 Log-log or semi-log graph paper or computer and software for ELISA data analysis. A four-parameter logistic (4-PL) curve-fit is recommended. Proteintech data analysis website, <https://www.ptgcn.com/products/elisa-kits/>.
- 3.8 Microplate thermostatic shaker.

4. Kit Components and Storage

Microplate - 96 well microplate precoated an anti-tag antibody (8 well × 12 strips)	1 plate	Unopened Kit: Store at 2-8°C for 6 months or -20°C for 12 months. Opened Kit: All reagents stored at 2-8°C for 7 days. Please use a new standard for each assay.
Protein standard - 8000 pg/bottle; lyophilized	2 bottles	
Capture antibody (100×) - 60 µL/vial*	1 vial	
Detection antibody, HRP-conjugated (100×) - 60 µL/vial*	1 vial	
Sample Diluent PT 4B1 - 30 mL/bottle	1 bottle	
Detection Diluent - 15 mL/bottle	1 bottle	
Wash Buffer Concentrate (20×) - 30 mL/bottle	1 bottle	
Tetramethylbenzidine Substrate (TMB) - 12 mL/bottle	1 bottle	
Stop Solution - 12 mL/bottle	1 bottle	
Plate Cover Seals	4 pieces	

* Centrifugation immediately before use

5. Safety Notes

- 5.1 Avoid any skin and eye contact with Stop Solution and TMB. In case of contact, wash thoroughly with water.
- 5.2 Do not use the kit after the expiration date.
- 5.3 Do not mix or substitute reagents or materials from other kit lots or other sources.
- 5.4 Be sure to wear protective equipment such as gloves, masks and goggles during the experiment.
- 5.5 When using an automated plate washer, adding a 30 second soak period following the addition of Wash Buffer to improve assay precision

6. Sample Collection and Storage

- 6.1 Serum: Allow blood samples to clot for 30 minutes, followed by centrifugation for 15 minutes at 1000×g. Clear serum can be assayed immediately or aliquoted and stored at -20°C. Avoid repeated freeze-thaw cycles.
- 6.2 Plasma: Use EDTA, heparin, or citrate as an anticoagulant for plasma collection. Centrifuge for 15 minutes at 1000×g within 30 minutes of collection. The plasma can be assayed immediately or aliquoted and stored at -20°C. Avoid repeated freeze-thaw cycles.
- 6.3 Cell Culture Supernatant: Remove particulates by centrifugation for 5 minutes at 500×g and assay immediately or aliquot and store samples at ≤ -20°C. Avoid repeated freeze-thaw cycles.
- 6.4 Tissue Homogenates: The preparation of tissue homogenates will vary depending upon tissue type. Rinse tissue with 1 x PBS to remove excess blood, Then add 1 mL of 1 x PBS per 100 mg tissue, homogenized and stored overnight at -20°C. After two freeze-thaw cycles were performed to break the cell membranes, the homogenates were centrifuged for 5 minutes at 10000 x g at 4°C. Collect the supernatant, assay immediately or aliquot and store at -20°C. Avoid repeated freeze-thaw cycles.

7. Regent Preparation

7.1 Wash Buffer (1X): If crystals have formed in the concentrate, warm to room temperature and mix gently until the crystals have completely dissolved. Add 30 mL of Wash Buffer Concentrate(20X) to 570 mL deionized or distilled water to prepare 1X Wash Buffer.


7.2 Antibody Cocktail (1X): Dilute 100X capture antibody and 100X HRP-conjugated detection antibody using Detection Diluent prior to assay. Suggested 1:100 dilution: 50 μ L 100X capture antibody + 50 μ L 100X Detection Antibody, HRP-conjugated + 4,900 μ L Detection Diluent. Mix gently but thoroughly.

7.3 Sample Dilution: Different samples should be diluted with corresponding Sample Diluent, samples may require further dilution if the readout values are higher than the highest standard OD reading. Variations in sample collection, processing and storage may affect the results of the measurement.

Recommended Dilution for different sample types: 1:16 or 1:32 is recommended for mouse serum and plasma; 1:160 or 1:320 is recommended for rat serum and plasma; 1:2 or 1:4 is recommended for cell culture supernatant and mouse tissue homogenate.

7.4 Standard Serial Dilution:

Add 2 mL Sample Diluent PT 4B1 in protein standard.



Add # μ L of Standard diluted in the previous step	—	500 μ L	500 μ L	500 μ L	500 μ L	500 μ L	500 μ L
# μ L of Sample Diluent PT 4B1	2000 μ L	500 μ L	500 μ L	500 μ L	500 μ L	500 μ L	500 μ L
	"sd7"	"sd6"	"sd5"	"sd4"	"sd3"	"sd2"	"sd1"

8. Assay Procedure Summary

Bring all reagents to room temperature before use (Detection antibody, HRP-conjugated antibody can be used immediately). To avoid cross-contamination, change pipette tips between additions of each standard level, between sample additions, and between reagent additions. Also, use separate reservoirs for each reagent.

8.1 Preset the layout of the microplate, including control group, standard group and sample group, take out the required number of microplate strips and return excess strips to the foil pouch containing the drying reagent pack and reseal; store at 4°C immediately. Microplate strips should be used in one week.

8.2 Add 50 µL standard or sample to appropriate wells. To avoid high background always add samples or standards to the well before the addition of antibody cocktail.

8.3 Add 50 µL 1× Antibody Cocktail solution (refer to Reagent Preparation 7.2) to each well. Seal plate with cover seal and incubate at 37°C on a microplate thermostatic shaker set at 400 rpm for 1 hour (incubate at 37°C for 2 hours is recommended if thermostatic shaker is not available) .

8.4 Wash

1) Gently remove the cover seal. Discard the liquid from wells by aspirating or decanting. Remove any residual solution by tapping the plate a few times on fresh paper towels.

2) Wash 4 times with 1× Wash Buffer, using at least 350-400 µL per well. Following the last wash, firmly tap plates on fresh towels 10 times to remove residual Wash Buffer. Avoid getting any towel fibers in the wells or wells drying out completely.

8.5 Add 100 µL TMB substrate solution to each well, protected from light. Incubate at 37°C on a microplate thermostatic shaker set at 400 rpm for 15 to 20 minutes. (Substrate Solution should remain colorless until added to the plate.)

8.6 Add 100 µL Stop Solution to each well in the same order as addition of the TMB substrate. Note: Avoid skin and eye contact with the Stop solution.

8.7 Read results immediately on a microplate reader at a wavelength of 450 nm. If possible, perform a double wavelength readout (450 nm and 630 nm).

8.8 Data analysis: Calculate the average of the duplicate readings (OD value) for each standard and sample, and subtract the average of the zero standard absorbance. Construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis, four-parameter logistic curve-fit (4-PL) analysis is recommended. If the samples have been diluted, the fitting result must be multiplied by the dilution factor used.


Procedure summary



9. Validation Data

9.1 Standard curve

These standard curves are provided for demonstration only. A standard curve should be generated for each set of samples assayed.



(pg/mL)	O.D	Average	Corrected
0	0.0412 0.0411	0.04115	-
62.5	0.1258 0.1254	0.1256	0.08445
125	0.2112 0.2047	0.20795	0.1668
250	0.3492 0.3419	0.34555	0.3044
500	0.5854 0.5634	0.5744	0.53325
1000	1.0271 1.0174	1.02225	0.9811
2000	1.7519 1.6787	1.7153	1.67415
4000	2.418 2.4264	2.4222	2.38105

9.2 Precision

Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested 8 times on one plate to assess intra-assay precision.

Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in 16 separate assays to assess inter-assay precision.

Intra-assay Precision					Inter-assay Precision				
Sample	n	Mean (pg/mL)	SD	CV%	Sample	n	Mean (pg/mL)	SD	CV%
1	8	2,032.5	38.5	1.9	1	16	2,074.0	103.6	5.0
2	8	518.4	13.2	2.5	2	16	518.2	9.7	1.9
3	8	257.8	5.5	2.1	3	16	256.1	7.0	2.7

9.3 Recovery

The recovery of mouse/rat HGF spiked to three different levels throughout the range of the assay in various matrices was evaluated.

Sample Type		Average% of Expected	Range (%)
Mouse serum	1:128	124	122-128
Rat serum	1:640	123	119-126
	1:1280	119	115-121
Cell culture supernatant	1:2	90	88-93
	1:4	123	118-127
Mouse tissue homogenate	1:8	118	116-119
	1:16	123	119-128

9.4 Sample values

Mouse/Rat serum - Mouse/Rat serum samples were evaluated for the presence of HGF in this assay.

Sample Type	Mean (ng/mL)	Range (ng/mL)
Mouse serum (n=16)	10.7	5.2-16.0
Rat serum (n=16)	97.6	53.1-238.7

Cell culture supernatant:

3T3-L1 cells were cultured in DMEM supplemented with 10% fetal bovine serum, 2.5 mM L-glutamine, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. An aliquot of the cell culture supernatant was removed, assayed for levels of HGF.

Mouse liver tissue, mouse kidney tissue and mouse spleen tissue from one adult mouse (chopped in 1-2 mm pieces in 30 mL of medium) was cultured for 3 days in RPMI supplemented with 10% fetal calf serum and 50 µM β-mercaptoethanol. Aliquot of the cell culture supernatant was removed, assayed for levels of HGF.

Rat spleen tissue from one adult rat (chopped in 1-2 mm pieces in 30 mL of medium) was cultured for 2 days in RPMI supplemented with 10% fetal calf serum and 50 µM β-mercaptoethanol. Aliquot of the cell culture supernatant was removed, assayed for levels of HGF.

	HGF (pg/mL)
3T3-L1 cell culture supernatant	2,968.8
Mouse liver cell culture supernatant	375.3
Mouse kidney cell culture supernatant	168.4
Mouse spleen cell culture supernatant	120.6
Rat spleen cell culture supernatant	131.0

Mouse tissue homogenate - Mouse kidney tissue, mouse spleen tissue and mouse liver tissue was rinsed with PBS to remove excess blood, chopped into 1-2 mm pieces, homogenized in 5-10 mL of PBS in a tissue homogenizer, and stored at ≤ -80°C 5 min. After two freeze-thaw cycles were performed to break the cell membranes, the homogenates were centrifuged for 5 minutes at 5000 x g to remove particulate. Homogenate from mouse kidney, spleen and liver were assayed for levels of mouse HGF.

	Mouse HGF (pg/mL)
Mouse liver tissue homogenate	3,366.5
Mouse kidney tissue homogenate	2,088.9
Mouse spleen tissue homogenate	2,070.6

9.5 Sensitivity

The minimum detectable dose of mouse/rat HGF is 7.8 pg/mL. This was determined by adding two standard deviations to the concentration corresponding to the mean O.D. of 20 zero standard replicates.

9.6 Linearity

To assess the linearity of the assay, samples were diluted with the appropriate **Sample Diluent** to produce samples with values within the dynamic range of the assay.

(The mouse serum was initially diluted 1:8. The rat serum was initially diluted 1:80.)

		Mouse serum	Rat serum	Cell culture supernatant	Mouse tissue homogenate
1:2	Average% of Expected	100	100	100	100
	Range (%)	-	-	-	-
1:4	Average% of Expected	102	94	95	102
	Range (%)	96-108	93-94	95-96	93-111
1:8	Average% of Expected	94	89	101	105
	Range (%)	94-95	85-93	99-104	89-121
1:16	Average% of Expected	84	90	97	104
	Range (%)	83-85	89-91	93-102	87-121

9.7 Specificity

This assay recognizes natural and recombinant mouse/rat HGF.

The following factors prepared at 50 ng/mL were assayed and exhibited no cross-reactivity or interference.

Recombinant mouse:

EGF

10. References

1. Szpirer C.et al. (1992). Genomics.13(2):293-300.
2. Seidel C.et al. (1992). Blood.91:806-812.
3. Oda Y.et al. (2000). HumPathol.31:185-192.
4. Margarida C.et al. (2003). Nat Med.9(11):1363-9.
5. Sheen-Chen S-M.et al. (2005). Cancer Epidemiol Biomarkers Prev.14:715-717.
6. Xie Q.et al. (2012). Proc Natl Acad Sci USA.109:570-575.
7. Conor AB.et al. (2017). Nat Rev Clin Oncol.14(9):562-576.
8. KimKH.et al. (2017). Exp Mol Med.e307(2017).
9. HaiyuWang.et al. (2020).Front Cell Dev Biol.8:55.