

For Research Use Only



# Atlantic Blue-conjugated Mouse IgG2a isotype control Monoclonal antibody

Catalog Number: PB-66360-2

## Basic Information

<b>Catalog Number:</b> PB-66360-2	<b>GenBank Accession Number:</b> GeneID (NCBI):	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 200 µg/ml by Nanodrop;	<b>Full Name:</b>	<b>CloneNo.:</b> 11A1B2
<b>Source:</b> Mouse		<b>Excitation/Emission maxima wavelengths:</b> 401 nm / 452 nm
<b>Isotype:</b> IgG2a		

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Megathura crenulata

## Background Information

This antibody is used as the isotype control of mouse IgG2a. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

## Storage

**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data