

For Research Use Only

FcZero-rAb® PE-CoraLite® Plus 594 Anti-Human CD3 (OKT3) Rabbit Recombinant Antibody

Catalog Number: PCL594-FcA65569



Basic Information

Catalog Number:
PCL594-FcA65569

Size:
100 tests, 5 ul/test

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC049847

GeneID (NCBI):
916

ENSEMBL Gene ID:
ENSG00000198851

UNIPROT ID:
P07766

Full Name:
CD3e molecule, epsilon (CD3-TCR
complex)

Calculated MW:
207 aa, 23 kDa

Purification Method:
Protein A purification

CloneNo.:
OKT3

Recommended Dilutions:
FC: 5 ul per 10⁶ cells in a 100 µl
suspension

Excitation/Emission maxima
wavelengths:
496 nm, 565 nm / 604 nm

Applications

Tested Applications:
FC

Species Specificity:
human

Positive Controls:

FC : human PBMCs,

Background Information

CD3 is a multimeric protein associated with the T-cell receptor (TCR) to form a complex involved in antigen recognition and signal transduction (PMID: 15885124). CD3 is composed of CD3γ, δ, ε, and ζ chains (PMID: 1826255). It is expressed by thymocytes in a developmentally regulated manner, T cells, and some NK cells (PMID: 3289580). The TCR recognizes antigens bound to major histocompatibility complex (MHC) molecules. TCR-mediated peptide-MHC recognition is transmitted to the CD3 complex, leading to the intracellular signal transduction (PMID: 11985657).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with FcZero-rAb™ CoraLite® Plus 488 Anti-Human CD19 (4G7), and 5 ul PE-CL594 Plus Anti-Human CD3 (OKT3) Rabbit IgG RecAb (PCL594-FcA65569, Clone: OKT3) or FcZero-rAb™ PE-CoraLite® Plus 594 Rabbit IgG Isotype Control Recombinant Antibody (PCL594-FcA98136, Clone: 240953C9). Cells were not fixed. Lymphocytes were gated.



1x10⁶ human PBMCs were surface stained with 5 ul PE-CL594 Plus Anti-Human CD3 (OKT3) Rabbit IgG RecAb (PCL594-FcA65569, Clone: OKT3) (red) or FcZero-rAb™ PE-CoraLite® Plus 594 Rabbit IgG Isotype Control Recombinant Antibody (PCL594-FcA98136, Clone: 240953C9) (blue). Cells were not fixed. Lymphocytes were gated.