

For Research Use Only



FcZero-rAb® PE-CoraLite® Plus 594 Anti-Human CD64 (10.1) Rabbit Recombinant Antibody

Catalog Number: PCL594-FcA65572

Basic Information

Catalog Number:
PCL594-FcA65572

Size:
100 tests, 5 µl/test

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC032634

GeneID (NCBI):
2209

ENSEMBL Gene ID:
ENSG00000150337

UNIPROT ID:
P12314

Full Name:
Fc fragment of IgG, high affinity Ia,
receptor (CD64)

Calculated MW:
374 aa, 43 kDa

Purification Method:
Protein A purification

CloneNo.:
10.1

Recommended Dilutions:
FC: 5 µl per 10⁶ cells in a 100 µl
suspension

Excitation/Emission maxima
wavelengths:
496 nm, 565 nm / 604 nm

Applications

Tested Applications:
FC

Species Specificity:
human, monkey, non-human primates

Positive Controls:

FC : human PBMCs,

Background Information

Fcγ receptor comprise a multigene family of integral membrane glycoproteins that exhibit complex activation or inhibitory effects on cell functions after aggregation by complexed immunoglobulin G (IgG) (PMID: 17005690). CD64, also known as Fcγ R1A, is a high-affinity receptor for the Fc region of IgG. It is expressed by monocytes/macrophages, activated neutrophils, dendritic cells, and early myeloid cells (PMID: 23293080; 19642859; 7680917). CD64 functions in both innate and adaptive immune responses.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with CoraLite® Plus 488 Anti-Human CD11b (ICRF44), and 5 ul PE-CoraLite® Plus 594 Anti-Human CD64 (10.1) Rabbit IgG RecAb (PCL594-FcA65572, Clone: 10.1) or Rabbit IgG RecAb (PCL594-FcA98136, Clone: 240953C9). Cells were not fixed.



1x10⁶ human PBMCs were surface stained with 5 ul PE-CoraLite® Plus 594 Anti-Human CD64 (10.1) Rabbit IgG RecAb (PCL594-FcA65572, Clone: 10.1) (red) or Rabbit IgG Isotype Control RecAb (PCL594-FcA98136, Clone: 240953C9) (blue). Cells were not fixed.