

For Research Use Only

FcZero-rAb® PE-CoraLite® Plus 594 Anti-Mouse CD86 (GL1) Rabbit Recombinant Antibody

Catalog Number: PCL594-FcA65673



Basic Information

Catalog Number:
PCL594-FcA65673

Size:
100ug, 200 ug/ml

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC013807

GeneID (NCBI):
12524

UNIPROT ID:
P42082

Full Name:
CD86 antigen

Purification Method:
Protein A purification

CloneNo.:
GL1

Recommended Dilutions:
FC: 0.25 ug per 10⁶ cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:
496 nm, 565 nm / 604 nm

Applications

Tested Applications:
FC

Species Specificity:
mouse

Positive Controls:

FC : LPS treated mouse splenocytes,

Background Information

CD86 (also known as B7-2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and CTLA-4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with CTLA-4 negatively regulates T-cell activation and diminishes the immune response. (PMID: 7513726; 1847722; 11029388; 27591335)

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 untreated (left) or LPS-treated (right) mouse splenocytes were surface stained with 0.25 ug PE-CoraLite® Plus 594 Anti-Mouse CD86 (GL1) Rabbit RecAb (PCL594-FcA65673, Clone: GL1). Cells were not fixed.



1×10^6 LPS treated mouse splenocytes were surface stained with 0.25 ug PE-CoraLite® Plus 594 Anti-Mouse CD86 (GL1) Rabbit RecAb (PCL594-FcA65673, Clone: GL1) (red) or 0.25 ug Rabbit IgG Isotype Control RecAb (PCL594-FcA98136, Clone: 240953C9) (blue). 1×10^6 untreated mouse splenocytes were surface stained with 0.25 ug PE-CoraLite® Plus 594 Anti-Mouse CD86 (GL1) Rabbit RecAb (PCL594-FcA65673, Clone: GL1) (black, dashed). Cells were not fixed.

