

For Research Use Only

# PE Anti-Mouse CD86 (GL1)

Catalog Number: PE-65068

14 Publications



## Basic Information

Catalog Number:

PE-65068

Size:

100ug, 0.1 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC013807

GeneID (NCBI):

12524

UNIPROT ID:

P42082

Full Name:

CD86 antigen

Purification Method:

Affinity purification

CloneNo.:

GL1

Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Cited Applications:

IF, FC

Species Specificity:

mouse

Cited Species:

mouse

## Background Information

CD86 (also known as B7-2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and CTLA-4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with CTLA-4 negatively regulates T-cell activation and diminishes the immune response. (PMID: 7513726; 1847722; 11029388; 27591335)

## Notable Publications

Author	Pubmed ID	Journal	Application
Xiaojin Chen	36297436	Pharmaceutics	FC
Lian Zhu	36321642	J Mater Chem B	FC
Lin Zhang	39841217	Naunyn Schmiedebergs Arch Pharmacol	FC

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> untreated or LPS-treated mouse splenocytes were surface stained with 0.2 ug PE Anti-Mouse CD86 (PE-65068, Clone: GL1). Cells were not fixed.



1X10<sup>6</sup> untreated (dashed line) or LPS-treated mouse splenocytes were surface stained with 0.2 ug PE Anti-Mouse CD86 (PE-65068, Clone: GL1) (red), or 0.2 ug Isotype Control (blue). Cells were not fixed.